

CCT Checklist/Survival Guide

PRE-SHIFT RESIDENT/CCT TEAM CHECKLIST:

Trauma Bay

CCT Surgical/Trauma Instrument Trays (R-side shelves)
*Please make sure the proper number of kits are present.
These are time-sensitive instrument trays that need to be available in the trauma room at all times.*

- Cordis (4)
- Central line kits (6)
- Peritoneal Lavage (2)
- Tracheostomy (3)
- Thoracotomy (2)
- Chest Tube Tray (3)
- Chest Tubes (Check ALL 3 Size Compartments)
- Pleurevac Drain (3)
- OB Kits (3)

L-side shelves

- Hemostatic dressings
- Epistaxis Nasal Packing
- Reboa
- Pelvic Binders
- Arterial Lines

Items to be aware of that are in the Pyxis in CCT 1

- Surgical Cric kit and Size 6 and 8 DCT Trachs
- Tourniquets
- Medtronic Pacer Pulse Generator and Magnet
- Tono-Pen
- Ring Cutter
- Bougies, Size 3 and 4 Hyperangulated Blades, Rigid Stylets
- McGrath Videolaryngoscope with Disposable Traditional Mac 3 and 4 blades
- Stryker Needle
- Blakemore Tube
- EZ-IO Device, Needles

Glidescope

- Make sure it is plugged in and functioning
- Supplies: S3/S4 blade x 2 each, 2-4 rigid stylets, 4 Bougies
- Disposable Intubating BFlex Bronchoscopes (3) Behind Door to Medication/Pyxis Room in CCT

Orange Airway Box

- Make sure at least one per red medical code cart
- Can have extra one in CCT 1 and in medication room

Each bed should have:

- BVM, Suction, NC, OxyMask immediately available
 - o Do NOT open yankauer or tubing
 - o Extra sealed supplies are in the airway box

All CCT Rooms

Each bed should have:

- BVM immediately available
- Suction
 - Canister attached is the ONLY thing that should be set-up
 - Sealed yankauer/tubing available and obtained if needed
 - Turn on suction to make sure it works
- Nasal Cannula and OxyMask
- Extra sealed supplies are in the airway box

Ultrasound

- Make sure it is plugged in and charging
- Supplies: Gel, Tegaderm, Probe cover, Wipes.
NO IV SUPPLIES
- Please clean the ultrasound before and after each use.** No Ultrasound Gel should be left on the US Probes.

Important materials located in the Supermarket (across from CCT Trauma Bay Room 1):

Please familiarize yourself with the following procedural trays/kits. Be sure to note where they are in the supermarket in case they are needed during a CCT shift.

- Wayne Pneumothorax Trays
- Extra Pleurevac Drains
- Pericardiocentesis Trays
- Lumbar Puncture Kits
- Thoracentesis/Paracentesis Trays
- Transvenous Pacer Kit
- Extra Central Line/Cordis Kits
- MidLine Catheter Trays
- Double Lumen Hemodialysis Catheters (16/20/24 cm)
- Extra Arterial Lines
- Traction Splint

Important: Do NOT open or preload blades, tubes, or yankauer sets. If any equipment/supplies noted above is running low or is not stocked, please notify the ED Support Team, and let the Charge Nurse Know. For used surgical instrument trays and/or airway boxes, call ext. 3644, Central Sterile for replacements. Be sure to clean up after use of a supply/instrument tray, or airway box.