

CCT ROOMS

- Ensure each bed in Rooms 1-5 has a working Bag Valve Mask
- Ensure each bed in Rooms 1-5 has a bougie taped above it.
- Ensure each bed in Rooms 1-5 has working suction (canister to tube to yankauer).
- Ensure Level 1 infusers are ready to be used.
- Ensure Ultrasounds are plugged in near trauma bay.
- Ensure Glidescope is loaded and functioning near trauma bay.
- Ensure you have the equipment needed to quickly draw a full set of labs and put in an IV.

TRAUMA BAY 1 – RED AIRWAY CARTS

- Ensure top shelf contains size 6 to 8 ET tubes (hooked up with working lights) with 2 stylets.
- Ensure both airway carts have working BVMs, suction, bougies.
- Ensure the drawers below the top drawer are reasonable (nursing will stock, check though).

BLUE / WHITE AIRWAY CARTS

- They are mobile with one typically placed near Rooms 2-3 and another placed near Room 5.
- These carts are to be organized, stocked, and maintained **by the residents.**
- There are three (3) drawers on these carts arranged in the following manner:

Drawer 1:

- Scalpels (11 blades)
- Tongue blades
- ETCO2 detectors
- Silk tape
- Laryngoscope handles
- Miller 3 and 4 blades
- Mac 3 and 4 blades
- One pre-prepared 7.0 ETT with stylet
- One pre-prepared 7.5 ETT with stylet

Drawer 2:

- Nasal cannulas
- Non-rebreather masks
- LMA (size 3.5 and 4.5)
- Nasopharyngeal airways
- Lubricating jelly
- Oropharyngeal airways
- Betadine bottle
- Benzoin bottle

Drawer 3:

- ETT fastener devices
- Extra stylets
- 20cc syringes
- Extra unopened endotracheal tubes.

If you are missing any supplies in the airway cart go to the CCT Charge Nurse

There are no Shileys or Cric kits in blue carts. Both are in bottom drawer of Trauma Bay carts.