

## CCT ROOMS

Each bed should have:

- ☐ BVM
- ☐ Working suction (canister, tubing, yankauer), Turn it on to make sure it works
- ☐ NC and Oxi-Mask. Turn on oxygen source to make sure it works

## TRAUMA BAY

- ☐ Ultrasounds plugged in
- ☐ Glidescope - turn it on
  - ☐ S3 blade
  - ☐ S4 blade
  - ☐ Rigid stylet
- ☐ Equipment needed to draw a full set of labs and put in an IV
- ☐ Bougies are in Pyxis, you will need to order them in order to get one
- ☐ Learn how to set up Level 1 infusers. These are not being prepared beforehand but will be set up when indicated.

## RED AIRWAY CARTS

- ☐ Ensure both airway carts have BVM, working O2 & suction
- ☐ Top shelf:
  - ☐ ET tubes: size 6-9
  - ☐ 2 stylets
  - ☐ 2 handles (check them to make sure light works!)
  - ☐ 2 Mac blades (3,4)
  - ☐ 2 Miller blades
  - ☐ ETCO2 detector
  - ☐ 10 CC syringe
  - ☐ Lubricant
- ☐ Ensure the drawers below the top drawer are reasonable (nursing will stock, check though).

- ☐ **BLUE / WHITE AIRWAY CARTS** These carts are to be organized, stocked, and maintained **by the residents.**

### ☐ **Drawer 1:**

- ☐ Scalpels (11 blades)
- ☐ Tongue blades
- ☐ ETCO2 detectors
- ☐ Silk tape
- ☐ Laryngoscope handles
- ☐ Miller 3 and 4 blades
- ☐ Mac 3 and 4 blades
- ☐ ET tubes: sizes 6-9

### ☐ **Drawer 2:**

- ☐ Nasal cannulas
- ☐ Non-rebreather masks
- ☐ LMA (size 3.5 and 4.5)
- ☐ Nasopharyngeal airways
- ☐ Lubricating jelly
- ☐ Oropharyngeal airways
- ☐ Betadine bottle
- ☐ Benzoin bottle

### ☐ **Drawer 3:**

- ☐ ETT fastener devices
- ☐ Extra stylets
- ☐ 20cc syringes
- ☐ Extra unopened endotracheal tubes.

\*\*\*If you are missing any supplies in the airway cart go to the CCT Charge Nurse\*\*\*