

# CCT Checklist/Survival Guide

## PRE-SHIFT RESIDENT CHECKLIST:

### Trauma Bay

CCT Surgical/Trauma Instrument Trays (R-side shelves)  
*Please make sure the proper number of kits are present.  
These are time-sensitive instrument trays that need to be available in the trauma room at all times.*

- ☐ Cordis (4)
- ☐ Central line kits (6)
- ☐ Peritoneal Lavage (2)
- ☐ Tracheostomy (3)
- ☐ Thoracotomy (2)
- ☐ Chest Tube Tray (3)
- ☐ Chest Tubes (Check ALL 3 Compartments)
- ☐ Pleurevac Drain (3)
- ☐ OB Kits (3)

**If there are any kits missing or if any are used during your shift, let the charge nurse know so they can be replaced (call ext 3644, Central Sterile for replacing surgical trays).**

### L-side shelves

- ☐ Hemostatic dressings
- ☐ Reboa
- ☐ Pelvic Binders
- ☐ Arterial Lines

### Items to be aware of that are in the Pyxis in CCT 1

- ☐ Cric kit and Size 6 and 8 DCT Trachs
- ☐ Tourniquets
- ☐ TV Pacemaker, Medtronic Pulse Generator and Magnet
- ☐ Tono-Pen
- ☐ Ring Cutter
- ☐ Extra bougies and Videolaryngoscope Blades
- ☐ Rigid stylets
- ☐ Stryker Needle
- ☐ Blakemore Tube
- ☐ EZ-IO Device, Needles

### Glidescope

- ☐ Make sure it is turned on and plugged in
  - ☐ Supplies: S3/S4 blade x 2 each, 2-4 rigid stylets, 4 Bougies
  - ☐ Disposable Bronchoscopes (3) Behind Door to Medication/Pyxis Room
- Do NOT open and preload blades or tubes.**

### Orange Airway Box

- ☐ Make sure at least one per red medical code cart
- ☐ Can have extra one in CCT 1 and in medication room

### Each bed should have:

- ☐ BVM, Suction, NC, OxyMask
  - o Do NOT open yankauer or tubing
  - o Extra sealed supplies are in the airway box

### All CCT Rooms

Each bed should have:

- ☐ BVM available
- ☐ Suction
  - Canister attached is the ONLY thing that should be set-up
  - Sealed yankauer/tubing available and obtained if needed
  - Turn on suction to make sure it works
- ☐ Nasal Cannula and OxyMask

### Ultrasound

- ☐ Make sure it is plugged in and charging
  - ☐ Supplies: Gel, Tegaderm, Probe cover, Wipes.
- NO IV SUPPLIES**
- ☐ **Please clean the ultrasound before and after each use.** No Ultrasound Gel should be left on the US Probes.

### Important materials located in the Supermarket (across from CCT Trauma Bay Room 1):

Please familiarize yourself with the following procedural trays/kits. Be sure to note where they are in the supermarket.

- ☐ Wayne Pneumothorax Trays
- ☐ Pericardiocentesis Trays
- ☐ Lumbar Puncture Kits
- ☐ Thoracentesis Trays
- ☐ Transvenous Pacer Kit
- ☐ Extra Central Line Kits
- ☐ MidLine Catheter Trays
- ☐ Double Lumen Hemodialysis Catheters (16/20/24 cm)
- ☐ Arterial Lines

**\* if any equipment or supplies noted above is running low or is not stocked, let the Charge Nurse know ASAP! \***